


TURP EQUIPMENT / RESECTOSCOPE

TURP Working Element

- Working Element enables the resection of the prostate with an electrical activated wire loop .
- Working Element also called TURP Working Element, Resectoscope Working Element.



Features :







- Designed to provide excellent utility.
- Use with 4mm 30 Degree Telescope.
- Compatible with Rotatable Sheath.
- Accurate dimensions.
- Non toxic.
- Advanced configuration.

Types

- Active Working Element
- Passive Working Element

TURP EQUIPMENT / RESECTOSCOPE

TURP Working Element



Types of TURP Loops

Cutting Loops : Remove prostate tissue in the form of small chips.

Cold Knife : Make incisions of ureteral stenosis for upper urinary tract obstructions.

Collin's Knife : Make straight incisions in a very hard prostate tissue.

Roller Electrode : Vaporize the prostate tissue using high power.

Ball Electrode 3mm : Coagulate and smoothen the freshly resected surface of the prostate gland.

Ball Electrode 5mm : Used when the surface to be coagulated is large.

Features

- Cutting and Coagulation.
- Precise design.
- Anti-corrosive.
- Compact design.
- Smooth finish .

Compatible With

Karl Storz.
Richard Wolf.
Olympus.

Types

Single Steam
Double Steam